

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

A.C. Module



A.C. Module

It is useful in Front of Air Conditioner to provide a sterile, bacteria free & dust free Air of 0.3 microns in a particular temperature . It works on the Laminar Air Flow Principle .

- HEPA of 0.3 μ of an efficiency of 99.97% down to 0.3 μ media is of glass fiber & made of pleates back back & forth in an anodized aluminum frame.
- Pre filter of dry fiber washable type with the frame of on all the sides.
- Statically Balanced motor blower Assembly (Heavy Duty)
- Velocity at the output of HEPA is 90 + 20 FPM.
- Wooden / M.S Duly Power Coated / Stainless Steel 304 Body

Working Area Size :
2' X 2'
3' X 2'
4' X 2'.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in