

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Capsule filling machine (laboratory size)



A versatile lever operated semiautomatic machine for various operations with levers is order of

1. plates for putting of capsule by hand
2. Removal of caps
3. filling of powder by dispensing of brush pressing
4. fixing of caps and sealing
5. Injection of capsules.

Contact parts of stainless for 100 capsules of size 0 to 2 or 5 (anyone size in order).

-do- capacity of 300 capsules at time.

Capsule filling machine (laboratory size)

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in