

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Centrifuge with heating arrangement



Centrifuge with heating arrangement

Features:-

- Electrically Operated Centrifuge
- Speed 1400 RPM (approx)
- Capacity 8 Tests.
- Inbuilt Heater.
- Digital RPM meter.
- Digital Temperature Controller.
- Digital Timer.
- Push Brake.
- Interlocking safety feature.
- Transparent Acrylic Lid