

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Dessiccator



Dessiccator

Plain Desiccators are physically the same as Vacuum Desiccator. Only difference is the absence of stop cock.

Grease should be applied on the flanges to make it air tight.

Material : Polycarbonate / Polycarbonate

SIZE (mm)

150
200
250
300

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in