

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## DIGITAL MELTING POINT APPARATUS WITH SILICON OIL BATH



This is a compact instrument designed for accurate determinations of Melting Points of solid samples.

### SALIENT FEATURES:

- \* Cylindrical Silicon Oil Bath
- \* Built-in Magnetic Stirrer with electronic speed controller.
- \* Electronic Controller for adjusting the heating rate.
- \* Digital Display of Temperature.
- \* Glare-free back-ground light with adjustable light Intensity.
- \* Provision for holding the melting point display

### DIGITAL MELTING POINT APPARATUS WITH SILICON OIL BATH

#### Specifications:

##### Working size:-

Temperature Range	2° above Ambient - 275°C
Readability	1°C
Display	3 Digit LED
Heating Rate Control	Electronic
Stirrer Speed Control	Electronic
Temperature Sensor	PT-100

NOTE: OIL BATH FOR TEMPERATURE MORE THAN 275 AND LESS THAN 350°C IS ALSO AVAILABLE ON REQUEST.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)