

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Digital pH and Conductivity Meter



Digital pH and Conductivity Meter

Specifications:

Display	3½ Digit LED
pH	
Range	0 - 14 pH
Resolution	0.01 pH
Accuracy	0.01 ± 1 Digit
Input Impedance	> 1012 ohms
Temp. Comp.	Manual: 0 - 100°C
Slope Correction	80 - 120%
CONDUCTIVITY	
Range	0 - 20 mS/Cm in 2 ranges
Resolution	1 S/Cm
Accuracy	± 0.5 % of Range ± 1 digit
Oscillator	1 KHz in-built
1 KHz in-built	0.2 - 1.8 adjustable directly on digital display
Temperature Compensation	0 - 50 °C adjustable from Temp. Table at 25 °C
TEMPERATURE	
Range	0 - 100 °C
Resolution	0.1 C
Accuracy	± 0.2% ± 1 digit
Sensor	Semi Conductor (PT-100) Type
Power Supply	220 V ± 10% AC, 50 Hz

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in