

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Digital Dissolution Rate Test Apparatus (Six Stage)



The Apparatus is same as the above but for six stages. The unit is equipped with, acrylic water bath controlled by digital temperature controller which is normally set at 37°C and duly fitted with pump which recirculates the water thus maintaining homogeneous temperature. The acrylic water bath is The digital timer can be set ranging from 1 to 999 minutes add it "BEEPS" at the end of the set time.

The unit is equipped with Digital setting of RPM as an optional which can be controlled from 25 to 150 RPM. The equipment is duly fitted with motorized lifting to facilitate easy operation of the stirrers' assembly. Digital dissolution rate test apparatus is equipped With

1. Digital RPM Indicator - Range 25 to 150 RPM
2. Digital Timer – Range 1 to 999 minutes
3. Digital temperature indicator – Range 30°C to 50°C.

Digital Dissolution Rate Test Apparatus (Six Stage)

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)