

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Distillation Apparatus (S.S.)-Table Pattern



- Double Walled Construction.
- Inner & Outer S.S. 304. The gap in between being filled with glass wool.
- The condenser is made of S.S. to work on 220/230 volts
- Delux Model Heavy Duty with Argon Welding

Capacity in Ltrs.

- 2 Litre/hrs
- 4 Litre/hrs
- 6 Litre/hrs
- 8 Litre/hrs
- 10 Litre/hrs

Optional:

Water Level Cut off System for safety of heaters in case of No / Insufficient water.

Distillation Apparatus (S.S.)-Table Pattern