

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Heating Block



Heating Block

This is a solid block heated unit with provision for different samples at a time. All the sample are heated at 150oC within $\pm 1\%$ accuracy with utmost uniformity of temperature These are fitted with Digital Temperature Controller cum Timer. Outer Body M.S. duly powder coated with Aluminium Block

Available Models:

- 6 Samples
- 15 Samples
- 20 Samples
- 25 Samples

OR AS PER REQUIREMENT WITH TUBE SAMPLES