

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Abel's Flash Point Apparatus



Abel's Flash Point Apparatus

Suitable for determination of close cup flash point determination of Petroleum and petroleum products having flash points below 70° C.

It comprises of test flame arrangement, oil cup, cover fitted with shutter mechanism, oil test jet or gas test jet, stirrer, water bath arrangement-temperature controlled by energy regulator and stand.

The unit is workable at 220 V AC supply.

OPTIONAL

Electrically Driver Stirrer with a FHP Motor having flexible shaft.

Spare Oil Cup Thermometer (-35°C to 70°C)

Water Bath Thermometer (-30°C to 80°C)

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in