

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Fruit Pressure / Firmness Tester



Fruit Pressure / Firmness Tester



DIGITAL MODEL

Fruit Firmness Tester also known as a Fruit Pressure Tester or Penetrometer measures fruit firmness and provides a quick, easy and accurate method to determine fruit maturity.

### ANALOG MODEL

Range 0-20 Kg / 44 Lb

Sensitivity 0.5 Kg.

Accuracy ± 2 graduations through 2500 Grams

### DIGITAL MODEL

Capacity 20 Kg / 44 Lb. / 196.10 Newton

Accuracy 0.5% +2 digits

Min. display 0.01Kg/0.01Lb/ 0.05Newton

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)