

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Gel Rocker / Shaker



Gel Rocker / Shaker

Can be used for many applications like gel staining, culture aeration, bathing, blotting membranes etc.  
A metal platform with non-slip rubber mat moves in a gentle 'See Saw' motion.  
The movement is controlled precisely with a DC motor.  
Cushion and rubber feet prevent unit from 'slipping'.

Motor	:	D.C.
Maximum RPM	:	20
Speed controller	:	Analogue
Platform Size	:	25 x 25 cm, 30 x 30 cm
Input Voltage	:	230V 50Hz, A.C. Supply.

**Optional:**

Digital RPM Indicator