

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Hardness Tester (Tablet )



Tablet Hardness Tester (Phyzer Type)



Tablet Hardness Tester (Monsanto Type)

A handy, tough and direct dial reading type tablet hardness tester, working as pair of pliers, generates breaking force in single grip. The instrument offers accurate and reproducible result much faster other conventional hardness tester.

### SPECIFICATIONS :

ACCURACY : + 0.2 kg  
LEAST COUNT : 0.2 kg  
SAFE TESTING FORCE : 20 kg  
MAXIMUM TABLET DIAMETER : 25 mm

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)