

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Hot Plate



Hot Plate

Rectangular / Round Hot Plate with outer M.S. sheet body duly powder coated/ S.S. 304. To work on 220/230 volts A.C. supply useful for continuous heating up to 350°C.

Available in models

M.S. Body with cast iron top / M.S. Body with S.S. top / S.S. Body with S.S. top

Size Available:

- a) 8" Dia round
- b) 10" x 12" x 6" (250x300x150mm) Rectangular
- c) 16" x 10" x 6" (400x250x150mm) Rectangular
- d) 18" x 24" x 6" (450x600x150mm) Rectangular

Available with Temperature controlling device such as

Energy Regulator

Or

Digital Temperature Controller

Or

Microprocessor Based PID controller



Hot plate Rectangular M.S body



Hot plate S.S. Body

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in