

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## LEAK TEST APPARATUS



LEAK TEST APPARATUS

It's the most ideal 'LEAK TEST' instrument for checking Blister Packs and Strips. It makes the conventional dessicator with a huge vacuum pump, obsolete. The German make Vacuum Pump ensures that a high level of vacuum is created in the shortest possible time. And the die - pressed. Polycarbonate vacuum dessicator sustains the vacuum for a long time. A vacuum gauge is provided to indicate the vacuum level. This is connected to the isolation valve which helps disconnect the source, hereby avoiding the need for running the vacuum pump continuously. A 0 – 5 minute timer helps the operator to attend to other chores at a time.

### FEATURES:

- Extremely Compact Construction
- Long Lasting, German Make Vacuum Pump
- Polycarbonate Vacuum Dessicator
- Stainless Steel (GMP) Construction
- Vacuum Level Upto 22" of Hg

### SPECIFICATIONS:-

Dessicator Diameter	250 mm(10")	300 mm (12")
Rating	5A 230V Single Phase	5A 230V Single Phase
Vacuum Level	22" Hg	22" Hg
Timer	5 minute (mechanical/Digital)	5 minute (mechanical/Digital)
Dessicator	Polycarbonate	Polycarbonate
Body	Stainless Steel	Stainless Steel

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)