

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Microscopic Stage warmer (For I.V.F. Labs)



For maintaining the temperature of critical biological specimens at desired level during microscopic observation. An anodized aluminum plate of less than 10 mm thickness with 50 mm central groove with glass and 175 (L) X 155(W) mm dimensions, which fits on most microscope stages Heated at desired temperature by a separate digital temperature display - cum controller with 2 meter long cable for convenient location.

The temperature range is from ambient to 55°C with 0.2°C accuracy.

Microscopic Stage warmer (For I.V.F. Labs)

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in