

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Orbital Shaker



Orbital Shaker (Heating Only)

Castor wheel mounted cabinet is double walled with inner chamber of stainless steel and outer casing of M.S. sheet duly painted / S.S. The gap of 75mm between the outer and the inner wall is filled with special grade glass wool to prevent thermal losses. Back of the unit is fitted with air circulation fan for maintaining the temperature uniformly throughout the chamber. Heating elements of nichrome / SS wire are placed in the air path and are duly insulated from the body.

(A) INCUBATOR SHAKER (Orbital Shaker Heating Only)

With Variable speed and Digital Temperature Controller
Inner Stainless Steel & Outer M.S. / S.S. with Acrylic Hood
Temperature range ambient to 700C.

Specially designed to hold :

- (i) To Hold 9 flasks of 250 ml OR Mix
(Size of platform 450 x 450 mm)
- (ii) To Hold 12 flasks of 250 ml OR Mix
(Size of platform 600 x 600 mm)

Optional :

Digital RPM Indicator
Microprocessor Based PID Controller

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in