

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Rotap Sieve shaker



Rotap Sieve shaker

Heavy duty(Improved Model ) FOR 8 ",12" & 18" TEST SIEVES The Rotap Sieve shaker is durably constructed and is run by ¼ H.P. Electric MOTOR 220 Volts A.C through reduction gear. Designed to carry upto 6/7 TEST SIEVES of 200 mm Dia (8") with one set lid & pan. It reproduces circular and hammering motion. Revolution approx.: 280 to 320 R.P.M and tapping approx.: 140 R.P.M to 160 R.P.M No foundations are required.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)