

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Table Top (Front Loading) Autoclave



- Front Loading Autoclave is fully automatic front loading autoclave with dry cycle.
- Electronic Circuit inside the autoclave with programmable integrated circuit chip.
- With dry cycle after sterilization without opening the door
- Digital Display of all the operations in sequence and timer & Temperature.
- The internal chamber have two number of trays made of stainless steel
- One Sterilization Cycle with dry cycle completes with single touching of start button.
- Fitted with Silicon Gasket, Inner Aluminium & Outer M.S. duly powder coated

Table Top (Front Loading) Autoclave

Inner Chamber Size	Ltrs (Approx.)
225 x 310 mm	12
255 x 310 mm	16
310 x 310 mm	22

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)