

SINCE 1978



OPTICS TECHNOLOGY



Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Digital Turbidity Meter



Digital Turbidity Meter

| Range | Normal Model | Deluxe Model |
|---------------|--------------------------------------|-----------------------------------|
| | 0-1000 (in 2 Ranges) | 0 - 1000 NTU (in 4 Ranges) |
| | 0-200 & 0-1000 | 0 - 2, 20, 200, 1000 NTU |
| Resolution | 0.1 NTU in 0-200, | 0.01 NTU in 0-2 & 0-20 |
| | 1 NTU in 0-1000 | 0.1 NTU in 0-200, 1 NTU in 0-1000 |
| Accuracy | ± 3% FS, ± 1Digit | |
| Display | 3½ Digit 7-Segment LED | |
| Response Time | Less than 10 seconds | |
| Light Source | 6V, 0.3Amp Tungsten Lamp | |
| Detector | Photodiode | |
| Sample System | 30 mm Clear Glass Test Tubes | |
| Accessories | Test tube Set of 5, Operation Manual | |

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in