

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Vacuum Dessiccator



Vacuum Dessiccator

Desiccators are made of Polypropylene and Polycarbonate. These autoclavable vacuum desiccators can hold vacuum upto 740 mm of Hg for 24 hrs. without any greasing.

The top dome, moulded in rigid and transparent Polycarbonate, gives a crystal clear view of the desiccant placed inside.

The Knurled knob on the top provides easy handling of the dome.

The lower part is made of Polypropylene. An internal groove is provided on the flange to hold a silicon rubber ORing.

The Polypropylene stopcock is fitted with a self lubricating PTFE plug.

This plug works three way i.e. vacuum creation, shutting off & vacuum releasing.

### SIZE (mm)

150  
200  
250  
300

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)