

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Vacuum Oven (Cylindrical)



Vacuum Oven (Cylindrical)

Inside chamber made of S. S. (304)
Outside Mild Steel Duly powder coated / S.S. 304 with Toughened Glass Door.
It is provided with 1 or 2 No. of S.S. trays.
Temp. Range Ambient Temp. +5°C to 130°C.
Temp. Controlled by Thermostat.
Vacuum Pressure indication by vacuum gauge.

Working Size:-

Dia X Depth	(Ltrs. Approx.)
9" X 12"	12
12" X 12"	21
12" X 18"	32
12" X 24"	43
15" X 20"	55
15" X 22"	62

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in