

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

VDRL Shaker



VDRL Shaker

Platform Size: 300 x 300 mm
Variable speed regulator for shaking speed upto 180 RPM
Supplied with cord and plug, operates on 230 V AC 50 Hz supply

OPTIONAL

Mechanical Timer (0-60 min)
OR
Digital Timer

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in