

SINCE 1978



# OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED  
An ISO 10002:2004 CERTIFIED

**CE Certified**

## Vortex Mixer / Vortex Shaker



Vortex Mixer / Vortex Shaker

Useful for accelerated mixing of solutions in test tubes, small flasks or centrifugal tubes etc.. The vibrating speed is regulated by a speed control device. Short stroke motor imparts agitation to the tubes or flasks. The equipment is workable on 220 V AC 50 HZ single phase. It is available with and without touch plate.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in), [opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website: [www.opticstechnology.in](http://www.opticstechnology.in)