

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Water Bath (Circulating)



Water Bath (Circulating)

This water bath is used to circulate hot water in a particular area / instrument like B.R. Apparatus OR ABBE Refractometer.

Temperature is controlled by Thermostat.

Inner Chamber SS-304, Outer MS duly powder coated/S.S.

Size of water bath : 8" x 8" x 8".

Temperature range: Ambient +5°C to 90°C.

Provided with variable speed regulator to control the speed of circulation.

Optional :

Digital temperature controller with an accuracy of $\pm 1^{\circ}\text{C}$

OR

Micro Processor based P.I.D. Controller with an accuracy of $\pm 0.5^{\circ}\text{C}$.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in