

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : LABORATORY TESTING INSTRUMENTS



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Serological Water Bath



Serological Water Bath

Double walled outer made of M.S. sheet duly painted / S.S. Inner stainless steel 304 grade sheet. The bath consist of 2 pilot lamps & on-off switch. Provided with immersion heater for heating.

Temp. controlled by Thermostat for controlling temp. from ambient temp.+5°C to 95°C with accuracy of $\pm 2^\circ\text{C}$

WORKING SIZE:-
300X255X175 mm
300X300X200 mm
300X300X300 mm
450X300X200 mm
Outer Body of construction :- M.S. duly powder coated body OR Stainless steel S.S. 304 Body

OPTIONAL

Low Water level cutoff device

Stirrer

Digital temperature controller with an accuracy of $\pm 1^\circ\text{C}$

OR

Microprocessor based P.I.D. Controller with an accuracy of $\pm 0.5^\circ\text{C}$.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in