

SINCE 1978



OPTICS TECHNOLOGY

Manufacturer & Exporters of : **LABORATORY TESTING INSTRUMENTS**



An ISO 9001:2008 CERTIFIED
An ISO 10002:2004 CERTIFIED

CE Certified

Tablet Making Machine (Motorised)



Tablet Making Machine (Motorised)

TABLET MAKING MACHINE (MOTORIZED):

12" DIA, STAINLESS STEEL : The machine is worked by oil immersed gearbox with 1/4 H.P. Motor driven, or by Heavy duty DC Geared Motor. Air blower is operated independently and air is thermostatically controlled to give temperature 30Å° to 110Å°C. Switches for motor & blower are fixed on a panel with neon indicators. Also mains lead with 15 amp plug provided.

2/3, DSIDC, Community Works Centre, J-Block, Shakurpur, (Opp. Britannia Ind. Ltd.) Delhi-110034 (India)

Phone: 0091-11-27101383, 27184490, Fax: 0091-11-27108890

E-mail: info@opticstechnology.in, opticstechnology@gmail.com

Website: www.opticstechnology.in